



VIDEO VRA-2

Video VRA™ is the latest behavioral conditioning system compatible with audiological testing for the 6 to 60-month infant/child population. The system incorporates proprietary hardware architecture with state of the art visual media material and equipment including video DVD disc players, LCD displays, switchbox routing and associated media. This system is useful in hearing assessment of hard to test populations, such as toddlers and developmentally delayed individuals. The Video VRA adds flexibility beyond the traditional, mechanical VRA/COR systems since all reinforcement is derived from an endless source of nonrepeating disc material. The lack of video redundancy should maintain patient interest levels and yield vital responses when hearing testing is conducted with operant conditioning techniques. The disc material by nature is flexible enough to accommodate locally popular material or current media such as animated movies, cartoons and other, high visual impact and attention getting material. In addition, sound associated with the disc and monitor display can be included with the reinforcement. Furthermore, this sound can be varied in loudness and calibrated to a meaningful output level. Indexing to specific chapters or playback reference points is simple process of noting and proceeding to a displayed target number. In most cases, any point will suffice. The standard playback speed “real-time” can be varied for maintaining the visual interest. Slow motion, fast motion, stop/freeze motion and other variations may be included on demand. Remember that the media is digital and the playback variations are sharp and clear, much like a photograph. Due to the degree of available media material and the flexibility of DVD this system has a wider appeal to a greater age range. In fact, you can customize any aspect of your reinforcement system from store purchased disc to custom, “homemade”, picture-disc.

Components for Deluxe System with centering monitor.

- DVD player with LCD Monitor/control unit and Centering unit with LCD screen
- DVD disc (early childhood oriented disc)
- Control box for routing visual reinforcement to selected monitor (not pictured)
- Two 19-inch televisions/LCD monitors (left, right screens)
- Patch cords and articulating wall mounts (not pictured)

Host site requirements include electric outlet within test room. Installation charges are extra. Monitor size upgrade is available for extra cost. Pricing and Specifications are subject to change. October 20, 2010. REV. 010202010

**For more information: SPA Associates, 545 MetroPlace South #100, Dublin, OH 43017
(800) 886-3245 WEB: www.spaassociates.com**